

Electromyography

Electromyography is a seductive muse.
-- Carlo DeLuca, *J Appl Biomechanics*, 1997

EID424
Bioengineering Applications in Sports Medicine
Oct 11, 2021
Prof. Kremenec


1

Electromyography

- + EMG
- + Measurement of the electrical signals produced by muscle tissue
 - + Quantify muscle activity
 - + Very low-level signals
- + Applications
 - + Clinical
 - + Research

2

Many Factors Affect EMG



3


EMG Acquisition

- + Electrodes
- + Amplifiers
- + A/D conversion
 - + Typically 12 or 16-bit
- + Analysis
 - + Amplitude measures
 - + Timing measures
 - + Frequency measures

4

Olde School EMG

- + Often used polygraph equipment
- + Counting boxes on graph paper



5

EMG Amplifiers

- + Differential amplifier
- + Must be low-noise, high-gain
 - + EMG very small signal (~100's of μV)
- + Subject must be isolated from power supply!
 - + Battery powered
 - + Optical isolation
 - + Transformer
- + Expensive

6

EMG Electrodes

- + Typically Ag-AgCl alloy
- + 2 per muscle
 - + Differential amp
 - + + 1 ground for system [usually]
- + Surface
- + Invasive
 - + Needle
 - + Fine wire

7

Surface EMG Electrodes

- + Easiest to use
- + Non-invasive
 - + Does not require special training/certification
- + Samples entire muscle
 - + Picks up signals from motor units beneath electrodes
- + Relatively good reproducibility
- + Good skin prep essential!

8

Surface EMG Skin Prep

- + All things between electrodes and muscle distort signal
 - + Skin
 - + Fat
 - + Hair
 - + Oil
 - + Lotion
- + Minimize impedance between electrode and muscle

9

Surface EMG Skin Prep Procedure

- + Shave if necessary
- + Abrade skin
 - + Remove layer of dead skin (and some live)
 - + Exfoliate
- + Clean with alcohol pad
 - + Ouch
 - + Be nice and do this first
- + Apply electrodes

10

Surface EMG Electrodes

- + Many types/shapes/sizes
 - + Single
 - + Bar
 - + Array
- + Reusable
 - + Must be kept clean
 - + Require application of electrolyte gel to ensure good contact
 - + Tape/glue to skin

11

Surface EMG Electrodes

- + Disposable
 - + Use once
 - + Pre-gelled/glued
 - + More commonly used

12


Pre-amplified Electrodes

- + Surface electrode with amplifier built-in
 - + Active vs passive pickups for guitar
- + Amplify signal at source
 - + Bigger signal relative to noise picked up by electrode leads, artifact from wires moving
- + Requires power
 - + More cumbersome
- + Not cheap

13

Surface Electrode Orientation

- + Pair of electrodes must be placed along direction of muscle fibers
 - + Measure the AP as it propagates along the fibers
 - + Try to avoid innervation zone of muscle



14

Invasive EMG Electrodes

- + Not distorted by skin/fat/hair/sweat/etc.
- + Can pick up signals from deep muscles
- + Samples small number of motor units
 - + Single motor unit?
- + Requires certification
 - + Training/certification requires medical background

15

Invasive EMG Electrodes

- + Confirm correct placement by stimulating through electrode
 - + Remove and replace if not correct
 - + "Poke and hope" technique [DeLuca]
- + Needle electrodes
- + Fine-wire electrodes
- + Ouch

16

Needle Electrodes

- + Very thin needle
- + Typically "Teflon"-coated
 - + Slides in and out easily
- + Exposed conductor only at end
- + Mostly used for stimulation
 - + Used for acquisition only when little or no motion involved

17

Fine-Wire Electrodes

- + Very thin wires inserted into needle
- + Insulated except at tip
- + Tip bent, sticking out of needle
- + Needle inserted, removed
 - + Barbed end of wire remains

18

Pros and Cons

- | | |
|---|---|
| <ul style="list-style-type: none">+ Surface<ul style="list-style-type: none">+ Relatively reproducible+ Gives idea of function of entire muscle+ Easy to use+ Signal distorted by intervening tissue | <ul style="list-style-type: none">+ Fine-wire<ul style="list-style-type: none">+ Hard to duplicate placement+ Only samples a few motor units+ Requires training, long setup+ No intervening tissue+ People don't like holes in their skin |
|---|---|

19


EMG Testing Procedure

- + Pair of electrodes per muscle
 - + Differential amplifiers
 - + Orient along direction of muscle fibers
- + Stimulator?
 - + Apply strong enough current to cause nerve/muscle to depolarize, generating action potential

20

So What's It Look Like?

- + Surface EMG interference pattern




21

EMG Interference Pattern

- + Result of many action potentials
 - + Different amplitudes (intervening tissue)
 - + Out of phase
 - + Some add, some subtract
- + For signal processing purposes, can be treated as a stochastic (random) process
 - + Band-limited Gaussian noise
 - + Zero-mean, variance a function of amplitude
 - + More on this later


22

Origin of EMG Interference Pattern



23

Motor Unit Recruitment



24

Factors Affecting EMG

- + Skin prep
 - + Reduce skin impedance
 - + More signal to amplifier
- + Electrode distance
 - + Farther apart → larger amplitude
 - + Standardize!
- + Electrode size
 - + Larger electrodes → more motor units

25

Comparing EMG

- + Take great care
 - + Different muscles have different properties
 - + Fiber type composition
 - + Length/tension
 - + Force/velocity
 - + Etc...
 - + Side-side, day-day in same person difficult
 - + Differences in placement, orientation of electrodes
 - + Between people almost impossible
 - + As above
 - + Also differences in size, strength

26

EMG Normalization

- + Enables comparison
- + Normalize to maximal voluntary isometric contraction (MVC)
 - + Most common
 - + Easiest conceptually
 - + Activity should isolate muscle in question as best possible
 - + Really max effort?
 - + Any relation to task being performed?

27

EMG Normalization

- + Normalize to standard activity
 - + Gait
 - + Lift standard weight
- + Normalize to data being analyzed
 - + Max in that activity is 1.0

28


EMG Crosstalk

- + Pick up unwanted signal from nearby muscle
- + Difficult to identify
 - + Cross-correlation? Not really
- + Double differential amplifier
- + Minimize with proper electrode spacing, good skin prep

29

Double-Differential Technique

- + Crosstalk should be similar in each single-differential signal
- + Another difference should eliminate crosstalk



30

Clinical EMG

- + What neurologist/orthopaedist/physiatrist is referring to when referring a patient for an EMG
- + Evaluating something specific
 - + Nerve injury
 - + Muscle weakness
 - + Carpal tunnel syndrome

31

Clinical EMG

- + Usually involves testing of one muscle/set of muscles or nerve
- + Usually involves stimulation of a motor nerve

32

Nerve Conduction Studies

- + Evaluate health of nerve by measuring its conduction velocity
 - + Speed with which AP propagates along its axon
- + Typically:
 - + EMG electrodes placed on muscle
 - + Nerve innervating muscle stimulated at several locations

33

Nerve Conduction Studies

- + Stimulate motor nerve
- + Pick up response at muscle
 - + Compound muscle action potential (CMAP)
 - + M-wave


34

Nerve Conduction Studies

- + Measure distance between sites of stimulation
- + Based on distance between sites and latency of response, conduction velocity along nerve can be estimated
 - + Amplitude, duration also of interest
- + Must have at least two sites!
 - + Only one site means time spent for signal to cross neuromuscular junction and propagate across muscle included
 - + Much slower than nerve conduction speed

35

Nerve Conduction Studies



+ Median nerve

NORMAL VALUES (MOTOR)

Segment	Distal Latency or NCV	Amplitude (mV)	Reference
Wrist to abductor pollicis brevis	3.9 ± 0.37 (3.4–4.5) msec	11.8 (7–17)	[2, 3]
Elbow to wrist	49.0 ± 3.9 (45.1–54.4) m/sec	11.8 (7–17)	[2, 3]
Axilla to elbow	56.3 ± 5.1 (50.0–68.3) m/sec	11.8 (7–17)	[2, 3]
Erb's to axilla	65.1 ± 6.1 (57.1–76.2) m/sec	11.8 (7–17)	[2, 3]
Erb's to elbow	62.9 ± 6.0 (51.0–76.0) m/sec	11.8 (7–17)	[3, 4]

36

F-Wave

- + Stimulus travels distally to muscle and proximally to motor neuron cell body
 - + *Orthodromically* and *antidromically*, respectively
- + Orthodromic response is M-wave
- + Neuron can fire to activate muscle in response to antidromic stimulus
 - + Latency, amplitude, duration of interest

37

Hoffman-Reflex (H-Reflex)

- + Stimulate nerve
 - + Some axons being stimulated are sensory nerve fibers from *muscle spindles*
 - + Sense muscle fibers stretching
- + These stimulate motor neuron to activate muscle

38

H-Reflex

- + Measure of motor neuron, spinal cord excitability
- + Amplitude of reflex, latency useful
- + Commonly measured in soleus muscle in back of leg
 - + Long nerve makes it easier to measure
 - + In shorter nerve, can be swamped out by initial motor response


39

H-Reflex

- + Very low stimulus gives no response
- + Gradually, evoke H-reflex
- + At maximal H-reflex, start of M-wave response
- + As M-wave increases, H-reflex goes away


40

H-Reflex



41

H-Reflex Recruitment Curves



42


Research Applications of EMG

- + Measure activity of muscles during various activities
- + Relate amount/type/timing of activation to the activity being performed
- + Holy grail: relate measured EMG activity to force being exerted by muscle
 - + Pretty much impossible

43

Typical Surface EMG Tracing

- + Surface EMG interference pattern
- + More or less band-limited, Gaussian noise



44

EMG Measurements

- + Amplitude measurements
- + Timing measurements
- + Frequency measurements


45

EMG Amplitude

- + Larger amplitude means more motor units recruited (mostly)
- + Amplitude *qualitatively* related to force generated in muscle contraction
 - + Many other factors related to force generation
 - + Other motor units that are not measured
 - + Type of motor units (I, II)
 - + Other muscles acting synergistically
 - + Co-contraction of antagonist muscles
 - + Pennation angle of muscle

46

Pennation Angle



- + Angle of muscle fibers
 - + Greater cross-sectional area
 - + Force-generating capacity is proportional
 - + Exert greater force
- + Remember:
 - + Want electrodes placed along direction of muscle fibers in order to measure propagation of AP along fibers


47

EMG Amplitude

- + Raw EMG requires some form of processing
 - + Otherwise, average amplitude is 0
- + Typically
 - + Rectified
 - + Rectified and smoothed (RMS)
 - + Sometimes referred to as integrated EMG (iEMG)
- + Don't forget to normalize!!!!


48

EMG Amplitude: Rectification



49

EMG Amplitude: Not What Is Meant by Integrated EMG




50

EMG Amplitude: Smoothed, Rectified EMG

- + Rectify EMG
- + Smooth result by either
 - + RMS of every xxx pts
 - + Number of points should be related somehow to the activity being studied
 - + Low-pass filter
 - + Integrate (and reset integrator!)

51

EMG Amplitude: Soothed, Rectified EMG



52

EMG Timing Measurements

- + When is muscle activated?
 - + Does it turn on when it should?
 - + On entire time it should?
- + When are agonist/antagonist activated relative to each other?
- + Amount of smoothing *very* important
 - + Don't want to miss events, or trigger incorrectly
- + Volumes of papers written about best way to determine when muscle is turned on
 - + How to pick threshold?

53

EMG Frequency Content

- + Very different between
 - + Surface
 - + In-dwelling
- + Beware noise at 50/60 Hz
 - + Don't forget harmonics
 - + Lots of signal as well
 - + Many systems include notch filter

54

Wire/Needle EMG Frequency Content

- + Can range from 20 Hz-several kHz
- + Typically not much more than 1 kHz

55

Surface EMG Frequency Content

- + Typically 20-500 Hz
- + Almost entirely done by 300 Hz

56

EMG Filters

- + Highpass filter
 - + Typically 10-20 Hz
 - + No EMG frequency components below 10 Hz
- + Mainly to remove motion artifact
 - + Noise created by wires, electrodes moving

57

EMG Filters

- + Lowpass filter
 - + \leq half sampling frequency
 - + Ideally, well under
 - + Limited by frequency content of EMG signal
 - + Int Soc of Electrophysiology and Kinesiology (ISEK) recommendations:
 - + > 350 Hz for surface
 - + > 450 Hz for wire
 - + > 1500 Hz for needle
- + Anti-aliasing
- + Remove high frequency noise

58

EMG Frequency Measurements

- + Mean frequency
- + Median frequency
- + Derived from power spectral density (PSD) computed via fast Fourier Transform (FFT)

59

EMG Mean Frequency

- + Centroid of the first moment of the power spectrum:
 - + $\Sigma [f G(f)] / \Sigma [G(f)]$
- + Can be adversely affected by high-frequency noise


60

EMG Median Frequency

- + Frequency which divides area under power spectrum in half
- + Typically lower than mean frequency

61

EMG Mean and Median Frequency



62

Factors Affecting EMG Frequency Content

- | | |
|--|--------------|
| + Interesting stuff: | + Bad stuff: |
| + Conduction velocity | + Skin |
| + Muscle diameter | + Fat |
| + These are essentially the shape of the action potentials | + Electrodes |
| + Slight effect of firing rate | |

63

Beware: EMG is Not Stationary

- + If body is moving, muscle is moving under electrodes
 - + Source of signal changing
- + Over time, EMG changes with fatigue
- + Only short periods of isometric contractions may be considered stationary!
- + Can only legally apply Fourier Transform to stationary signals



64

Moving Screws Everything Up

- + With movement:
 - + Muscle fibers/motor units beneath electrodes move
 - + Change shape
 - + Change location relative to electrodes
- + At different joint angles more/less/different motor units may be active to maintain the same force level

65

EMG and Fatigue

- + Fatigue
 - + Inability to maintain the required force
- + Peripheral causes
 - + Metabolic processes
 - + Issues in muscle, neuromuscular junction
- + Central causes
 - + Failure of central nervous system to drive muscle
 - + Not motivation


66

EMG and Fatigue

- + In sub-maximal isometric contractions
 - + EMG amplitude *increases* with fatigue
 - + More motor units recruited to perform same work
 - + EMG mean/median frequency decreases with fatigue
 - + Decrease of AP conduction velocity along muscle fiber with build-up of metabolites
 - + Power spectrum compresses, becomes more skewed
 - + Motor unit synchronization

67

EMG and Fatigue



68


Fatigue

- + Women more fatigue-resistant than men
- + Elderly more fatigue-resistant than young
- + Differences in:
 - + Fiber type
 - + Strength
 - + Blood flow
 - + ...

69

Peripheral vs Central Fatigue


- + Central activation ratio (CAR)
- + Subjects performs maximal voluntary contraction, measure force
- + Supramaximal stimulus superimposed
 - + Additional force: muscle can be driven further; fatigue is failure of nervous system drive
 - + No additional force: failure in force generation is in the muscle



70

Cutting Edge [as in, still unexplained]*

- + Distance cycling
 - + Two hours
 - + Interspersed sprints (every 20 min), 3k time trial at end
- + Men: central and peripheral fatigue
- + Women: central fatigue only



*Glance, Kremenich, McHugh. Eur J Appl Physiol 2013 Apr;113(4):1091-8.

71

Strange EMG Changes...

	Men	Women	Sig?
VL amp	Unchanged	Unchanged	0.629
VM amp	Unchanged	↓ 31%	0.021
RF amp	Unchanged (↓ 12%)	↓ 17%	0.569
VL Mean freq	Unchanged	Unchanged	0.369
VM Mean freq	Unchanged	Unchanged	0.757
RF Mean freq	↑ 26%	Unchanged (↑ 5%)	0.097

- + Most research done with fatigue uses pretty high-level contractions
 - + Cycling is *extremely* low-level
 - + Different mechanism/motor unit rotation strategy

72


Newer EMG Application: Magnetic Stimulation

- + Magnetic stimulator
 - + Think Faraday's Law
 - + Induce current in wire with magnetic field around wire
 - + Replace "wire" with "nerve"
 - + Painless, especially compared with electrical stimulation
- + Stimulation of peripheral nerves
- + Stimulation of nerve roots at spinal level
- + Stimulation of motor cortex
 - + Recently, FDA-approved for treatment of depression

73

Transcranial Magnetic Stimulation

- + Stimulate motor cortex of brain
- + Measure response at muscle
- + Evaluate entire conduction path
- + Typically done using arm, hand muscles



74

Transcranial Magnetic Stimulation

- + Motor evoked potentials (MEP)
 - + Stimulate, measure size of AP response
 - + Measure *central motor conduction time* (CMCT)
 - + Also stimulate nerve root in spinal cord
 - + Often easier to evoke while subject is performing a contraction
 - + Cortex more *excitable*
 - + Site of fatigue?
- + Silent period after MEP during active contractions

75

Bathroom Reading

- + Basmajian, DeLuca CJ. *Muscles Alive*. Baltimore: Williams and Wilkins, 1985.
- + DeLuca CJ. The Use of Surface Electromyography in Biomechanics. *J Appl Biomech*, 1997, 13:135-163.
- + Kamen G, Caldwell GE. Physiology and Interpretation of the Electromyogram. *J Clin Neurophysiol*, 1996, 13(5):366-384.
